



Navtech Practice Solutions
1979 Stout Drive, Unit 2
Warminster, PA 18974
(215) 364-2012
www.navtech360.com



OWANDY-CAM^{HD}

SIMPLIFY COMMUNICATION WITH YOUR PATIENTS

USB DIGITAL INTRA-ORAL CAMERA

Owandy-Cam




Facilitates communication with the patient and acceptance of care

Explain clinical procedures to your patients more easily! Your patient feels more involved and more confident.



4 sockets of view

Its "focus free" will give you a great depth of field and high quality images. Requiring no adjustment, get a perfectly sharp image of the single tooth smile and save time!



SMILE

INTRAORAL

INTRAORAL

MACRO



100% integrated

Owandy-Cam HD is fully integrated into Julie and the main software on the market.

Owandy
RADIOLOGY

Owandy
RADIOLOGY



Form ergonomic and wide angle lens of 105 °

Its rounded shape and ultra-thin head improve patient comfort and prevent trauma in the mouth.

With a viewing angle of 105 °, it differs from other cameras on the market by allowing better exploration of distal areas.



Easy of use

Fix images with a simple swipe of the touch button. You won't have any blurry images and you can work with confidence.

Thanks to its USB 2 connection, you can easily move and share the camera with your colleagues and associates in all the rooms of your practice.

Characteristics techniques

Sensor	High sensitivity ¼ "Sony CCD
Resolution	50 frames / sec
Lighting	8 LEDs
Focus	Focus free large depth of field
Viewing angle	105 °
USB cable length	2.5 m
Handpiece size (L x W x H) Intraoral part	205 x 28 x 24 mm
thickness Handpiece weight	11 mm
	55 g
Image freeze	Touch button



~ Push & Pull system

Owandy-Cam HD is easily integrated in the chair like your handpieces.



Owandy-Cam HD

Owandy-Cam HD

DIGITAL WORKFLOW OWANDY RADIOLOGY

A COMPLETE RANGE TO MEET YOUR NEEDS



BEST Communication

